

Ready for a quick dip?

DRUGCHECK®

From the maker of the NxStep Onsite Test Cup

With our Dip Drug Test, you can now experience the same reliability and simplicity you've come to trust from the DrugCheck® name. We know dependability is your top priority. That's why we created an alternative test in a dip card format with the same accuracy and ease of use as the DrugCheck® NxStep Onsite Test Cup.

Easy, reliable, *fast*

The DrugCheck® Dip Drug Test is a rapid, qualitative immunoassay for the detection of up to 10 drugs, at the SAMHSA cutoff levels, in a convenient format. Simply insert the test strips in the collected urine sample for 10-15 seconds and read the results at 5 minutes.


- ✓ Accurate results in 5 minutes
- ✓ Easy to read
- ✓ Simple procedure
- ✓ 18-month shelf life
- ✓ Available in single- and multiple-test panels
- ✓ No urine to pipette
- ✓ No instrumentation
- ✓ Made in U.S.A.




Dip device into collected urine for 15 seconds



Replace cap



Read results at 5 minutes



AMP AMP	Amphetamine	1000 ng/ml
BAR BAR	Barbiturate	300 ng/ml
BZD BZD	Benzodiazepine	300 ng/ml
BUP BUP	Buprenorphine	10 ng/ml
COC COC	Cocaine	300 ng/ml
MET MET	Methamphetamine	1000 ng/ml
MTD MTD	Methadone	300 ng/ml
OPI OPI	Opiates	300/2000 ng/ml
OXY OXY	Oxycodone	100 ng/ml
PCP PCP	Phencyclidine (PCP)	25 ng/ml
THC THC	Marijuana	50 ng/ml
MDMA MDMA	Ecstasy	500 ng/ml
PPX PPX	Propoxyphene	300 ng/ml
TCA TCA	Tricyclic Antidepressants	1000 ng/ml

Preparation, Specimen Collection, Test Procedure, and Interpretation of Results

Collectors should be trained on the procedures for operating the tests. Do not break the seal of the protective pouch until ready to begin testing. Visually inspect the foil package to ensure it has not been compromised before beginning the test. If the package is not intact, the integrity of the test may be compromised.

1. Fresh urine specimens should be collected directly into the collection cup and do not require any special handling or pre-treatment.
2. The **DRUGCHECK®** Dip Test should be used with a collection cup that employs a temperature strip to validate that the urine is a freshly voided specimen. This device should be checked immediately after collection.
3. Require the donor to present a photo identification.
4. Ask donor to remove any unnecessary outer clothing and empty all pockets.
5. Keep all backpacks and purses out of the restroom. Insure that all backpacks, wallets and purses are secured, preferably in a locked cabinet. If using a lockable cabinet or container, make sure that the donor is given the key.
6. Secure the collection site restroom by doing the following;
 - a. Remove all cleaning supplies or substances that could be used to adulterate the specimen
 - b. Tape off any hanging soap dispenser

- c. Move waste can to the outside of the restroom
 - d. Secure water source (faucets and spouts) using masking or duct tape
 - e. Add bluing to tank and bowl water
7. Tear open the foil pouch and remove the Dip Test.
 8. Issue the collection cup to the individual to be tested. Instruct the donor not to flush the toilet or turn the water on and to open the door when they are done.
 9. Have the individual urinate directly into the collection cup. Ensure that the specimen contains at least 30 mL of urine.
 10. Check temperature (normal range between 90°-100°F or 32°-38°C) and observe specimen for foreign material and/or discoloration.
 11. Have available for use a Chain of Custody form (COC).
 12. Lower the Dip Test into the specimen, ensuring that the bottom of the test strips contact the urine. **Plastic dip device should not touch the specimen (urine).**
 13. Replace the cap covering the ends of the test strip.
 14. Read results at 5 minutes.

Drug Test Results

Presumptive Positive and Confirm: One colored line only is visible in each Control Zone. No colored line whatsoever appears in the appropriate Test Zone, indicating a presumptive positive result for the corresponding drug of that specific Test Zone. Send this urine specimen to a certified laboratory for confirmation.

Negative: Two colored lines are visible in each Control Zone and in the appropriate Test Zone, indicating that the concentration of the drug is below the detection limit of the test or the drug is not present.

Retest: No Lines - If a colored line is not visible in the Control Zone, the test is invalid. Another test should be run to re-evaluate the specimen. Each strip in the **DRUGCHECK®** Dip Test is read and functions independently. An invalid result on one test strip does not invalidate other results derived from the same device.

NOTE: There is no meaning attributed to line color intensity or width. Any evidence of a line should be considered a line. In order to prevent any incorrect results, the test results should not be interpreted after 10 minutes.


If confirmation is required, secure specimen in appropriate shipping container and send to Confirmation Lab with signed COC.



Dip device into collected urine for 15 seconds



Replace cap



Read results at 5 minutes